



Prusa Protective Face Shield - RC3.1 - stack versions



VIEW IN BROWSER

updated 31. 3. 2020 | published 27. 3. 2020

Summary

1h 20min per pcs

Stacking versions of Prusa Protective Face Shield - RC3

A prototype face shield that we developed. In three days, we went through dozens of prototypes and two verifications with the Czech Ministry of Health.

Provided G-codes are only for Original Prusa i3 MK3/S printers. Do not run the G-codes on other printers. Risk of damage.

We provided two options, including reinforcement parts:

- 4 pcs on the print bed (2 stacks, each with 2 headbands)
- 8 pcs on the print bed (2 stacks, each with 4 headbands)

Version RC3.1 changelog Changes from RC3

• Fixed missing supports under two pins in the second headband

The printer will be running at increased speeds. **Expect louder** operation and make sure your printer is in good shape before running this G-code. Check the Maintenance guide.











8





Healthcare > Accessibility

Supports are included in the model. Print with at least 2 perimeters, and 5% infill. The model is optimized for a layer height of 0.25mm.

Print it from PETG, headbands from other material can be brittle with these print settings.

This remix is based on



Prusa Protective Face Shield - RC3

by Prusa Research

Print Files (.gcode)

↓ DOWNLOAD ALL FILES

72×72

covid19_rc31_4pcs_complete_025mm_pe 17.7 MB

updated 27. 3. 2020

(5.31 hrs ≡ 0.25 mm ♣ 0.40 mm ♠ PET ₫ 179.00 g



72×72

covid19 rc31 8pcs complete 025mm pe 36.6 MB

updated 27. 3. 2020

updated 27. 3. 2020

① 11.11 hrs≡ 0.25 mm ♣ 0.40 mm � PET ଣିଁବି 361.10 g



Model Files (.stl, .3mf, .obj, .amf)

→ DOWNLOAD ALL FILES

72×72

covid19 rc31 4pcs complete.3mf

479.4 KB





License **G**





Attribution-NonCommercial

- **★** | Sharing without ATTRIBUTION
- ✓ | Remix Culture allowed
- **≭** | Commercial Use
- **≭** | Free Cultural Works
- **x** | Meets Open Definition